



EZ-fill® Cartridge Platform



Clean. Sterile. Ready to Be Filled.

Stevanato Group offers the broadest portfolio of EZ-fill® ready-to-use glass cartridges, designed to support the evolving needs of pharmaceutical manufacturers and advanced injectable therapies.

EZ fill® cartridges are washed, siliconized, capped and sterilized, enabling efficient filling operations while ensuring high quality, consistency, and regulatory compliance.

Available in multiple secondary packaging configurations to fit different fill & finish set-ups.

Designed for a Broad Range of Drug Delivery Systems

Engineered for seamless integration with pen injectors, autoinjectors, on-body delivery systems, and more,

EZ-fill® cartridges deliver reliable performance throughout the entire product lifecycle from filling to injection.



Features and Benefits.



NEXA® QUALITY

- Simplified and de-risked fill & finish operations by eliminating washing, siliconization and crimping steps
- No glass-to-glass contact throughout manufacturing and transport
- High cosmetic quality and mechanical resistance
- Reduced risk of rejects and breakages



REDUCE TCO

- Reduced capital investment through streamlined fill & finish operations
- Smaller manufacturing footprint and optimized workforce requirements
- Lower validation and maintenance burden
- Improved fill & finish yield by reducing drug product losses



INCREASE FLEXIBILITY


- Flexibility to leverage multi-purpose fill & finish lines for both cartridges and syringes
- Different primary containers can be processed on the same aseptic filling line without changing the infeed module




FASTER TIME TO MARKET

- Pre-validated solutions and pre-verified gliding profile to support optimal device performance
- Broad adoption across a network of strategic CDMO partners, enabling both development and commercial phases

Preferred Market Choice.




A PROVEN PROCESS



PACKAGING

Cartridges are delivered in nest and tub and tray configurations, in protective bag compatible with steam and EtO sterilization



The Right Choice for Every Fill & Finish Process.



PRE-CAPPED CARTRIDGE PORTFOLIO

FORMAT	OUTER Ø	INNER Ø	HEIGHT	FINISH	CLOSURES	PCS
1.5 mL	8.65	6.85	56.95	7.15 mm	Combi-seal	100
3 mL	11.6	9.7	62.3	7.15 mm	Combi-seal	100
5 mL long	14.0	12.0	65.0	7.15 mm	Combi-seal	80
5 mL	22.05	19.05	46.0	13 mm	Front Stopper & Alucap	42
10 mL	22.05	19.05	68.7	13 mm	Front Stopper & Alucap	42
20 mL	26.5	23.85	72.5	13 mm	Front Stopper & Alucap	36



UN-CAPPED CARTRIDGE PORTFOLIO

FORMAT	OUTER Ø	INNER Ø	HEIGHT	FINISH	CLOSURES	PCS
3 mL	11.6	9.7	62.3	7.15 mm	Un-capped	Nest & Tub 100
3 mL	11.6	9.7	62.3	7.15 mm	Un-capped	Tray 220

All formats are in Nexa® configuration.

The height refers to the glass cartridge without the assembled closure.

Headquarters

Via Molinella, 17
35017 Piombino Dese
Padova, Italy

**Global presence with operating units
and commercial offices in 10 countries.**

Visit our website to contact us.



stevanatogroup.com

SG, Stevanato, EZ-fill, EZ-fill Smart, Alba, Nexa, Fina, Alina, Vertiva are registered trademarks of Stevanato Group and/or its affiliates in one or more jurisdictions.

EZ-FILL_CAR_202609