



Ampoules

High-quality pharmaceutical glass containers with a broad range of geometries

Stevanato Group ampoules are a high-quality primary packaging solution for pharmaceutical applications. The glass ampoules are produced with Type 1 Glass under secured and certified process conditions in accordance with the current Pharmacopoeiae (EP, USP, JP), guaranteeing optimized technical and cosmetic performance.

Formats

To meet the different pharmaceutical needs, Stevanato Group production includes all formats of ampoules, classified as **"cut open"** and **"flame-cut open"** form B, **"open cup"** form C, **"stem sealed dome type"** form D and **"stem sealed Marzocchi type"** form E.

All formats are available with **color break**, **OPC** and **score ring**, their volumes range **from 1 ml to 25 ml**.



Form B



Form C



Form D



Form E
Marzocchi



Color
Break



Score
Ring



OPC




Ampoules

Benefits

Stevanato Group ampoules production methods and logistical processes are compliant with the **strict terms and conditions** of pharmaceutical and cosmetic markets in order to ensure the correct connection to the client's QA & QC processes.

QUALITY CONTROLS

Ampoules quality control covers every phase of the process. Production lines are provided with:



Equipment for the **breaking force measurement**, with automatic collection and analysis of data on host computers



Integrated system for data management of products quality control



Systems for the **traceability of all manufacturing parameters** (materials, equipment, operators, etc...)

LOGISTICS



Thanks to its automated logistics system, Stevanato Group ensures a punctual delivery of the product.

Options

Upon customer's request, Stevanato Group can also provide an **outer gliding treatment**, inner **ammonium sulfate treatment**, and **printing**.

