

Products & Capabilities

EZ-fill®
Cartridge Platform

Stevanato Group

MANAGING COMPLEXITY, DELIVERING VALUE

Stevanato Group

Clean. Sterile. Ready to Be Filled.

Stevanato Group offers the broadest portfolio of EZ-fill® ready-to-use glass cartridges, designed to support the evolving needs of pharmaceutical manufacturers and advanced injectable therapies.

EZ fill® cartridges are washed, siliconized, capped and sterilized, enabling efficient filling operations while ensuring high quality, consistency, and regulatory compliance.

Available in multiple secondary packaging configurations to fit different fill & finish set-ups.

Designed for a Broad Range of Drug Delivery Systems

Engineered for seamless integration with pen injectors, autoinjectors, on-body delivery systems, and more,

EZ-fill® cartridges deliver reliable performance throughout the entire product lifecycle from filling to injection.


Vertiva®
Alina®
Pen injector Platform
Deora™ 50

Features and Benefits.



NEXA® QUALITY

- Simplified and de-risked fill & finish operations by eliminating washing, siliconization and crimping steps
- No glass-to-glass contact throughout manufacturing and transport
- High cosmetic quality and mechanical resistance
- Reduced risk of rejects and breakages



REDUCE TCO

- Reduced capital investment through streamlined fill & finish operations
- Smaller manufacturing footprint and optimized workforce requirements
- Lower validation and maintenance burden
- Improved fill & finish yield by reducing drug product losses



INCREASE FLEXIBILITY


- Flexibility to leverage multi-purpose fill & finish lines for both cartridges and syringes
- Different primary containers can be processed on the same aseptic filling line without changing the infeed module




FASTER TIME TO MARKET

- Pre-validated solutions and pre-verified gliding profile to support optimal device performance
- Broad adoption across a network of strategic CDMO partners, enabling both development and commercial phases

Preferred Market Choice.




A PROVEN PROCESS



PACKAGING

Cartridges are delivered in nest and tub and tray configurations, in protective bag compatible with steam and EtO sterilization




Protective Bag
Available in single- or double-bag configurations, depending on the selected de-bagging process

Liner
Provides protection for nested containers during lid opening

Sealing Lid
Sealed to the tub and designed for an easy and smooth peel-off

Nest & Tub
Designed to support different filling technologies

The Right Choice for Every Fill & Finish Process.



PRE-CAPPED CARTRIDGE PORTFOLIO

FORMAT	OUTER Ø	INNER Ø	HEIGHT	FINISH	CLOSURES	PCS
1.5 mL	8.65	6.85	56.95	7.15 mm	Combi-seal	100
3 mL	11.6	9.7	62.3	7.15 mm	Combi-seal	100
5 mL long	14.0	12.0	65.0	7.15 mm	Combi-seal	80
5 mL	22.05	19.05	46.0	13 mm	Front Stopper & Alucap	42
10 mL	22.05	19.05	68.7	13 mm	Front Stopper & Alucap	42
20 mL	26.5	23.85	72.5	13 mm	Front Stopper & Alucap	36



UN-CAPPED CARTRIDGE PORTFOLIO

FORMAT	OUTER Ø	INNER Ø	HEIGHT	FINISH	CLOSURES	PCS
3 mL	11.6	9.7	62.3	7.15 mm	Un-capped	Nest & Tub 100
3 mL	11.6	9.7	62.3	7.15 mm	Un-capped	Tray 220

All formats are in Nexa® configuration.

The height refers to the glass cartridge without the assembled closure.

Headquarters
Via Molinella, 17
35017 Piombino Dese
Padova, Italy

Global presence with operating units and commercial offices in 10 countries.

Visit our website to contact us.

Stevanato Group

stevanatogroup.com

EZ-FILL_CAR_202609

, Stevanato, EZ-fill, EZ-fill Smart, Alba, Nexa, Fina, Alina, Vertiva are registered trademarks of Stevanato Group and/or its affiliates in one or more jurisdictions.